

Product Datasheet

CD154 antibody (PE) (orb308750)

Description

Mouse monoclonal antibody against CD154 conjugated to

Reactivity Human, Primate

Conjugation PE

Tested Applications FC

Immunogen human CD154 fusion protein

Target CD154

Preservatives Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note For research use only

Application notes Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype Mouse IgG1

Clonality Monoclonal


Purity Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID [P27548](#)


Entrez [21947](#)

Dilution Range Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date 12 months from date of receipt.



Separation of human CD154 positive lymph...



Flow cytometry surface staining pattern ...