


Product Datasheet

CSF3R antibody (orb308866)

Description	Rabbit polyclonal antibody to CSF3R
Species/Host	Goat
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	G-CSF receptor / CSF3R (aa551-563)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 92.1kDa according to NP_000751.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 1-3µg/ml.
Protein Sequence	PPELGKSPLTHYT
Clonality	Polyclonal
MW	92.1
Entrez	1441
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Peripheral Bloo...