


Product Datasheet

EPHX2 antibody (orb308868)

Description	Rabbit polyclonal antibody to EPHX2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EPHX2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in Human Kidney lysates (calculated MW of 57.1kDa according to NP_001243411.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	SRTFKSLFRASDES
Clonality	Polyclonal
MW	62.6; 57.1; 55.6
Entrez	2053
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...