


Product Datasheet

EXOG antibody (orb308871)

Description	Rabbit polyclonal antibody to EXOG
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	EXOG
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40kDa band observed in lysates of cell lines HepG2 and Jurkat (calculated MW of 41.1kDa according to NP_005098.2). Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 5µg/ml.
Protein Sequence	KAKEQSGTQIRKP
Clonality	Polyclonal
MW	41.1; 35.4
Entrez	9941
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 (Lane 1) ...



Immunohistochemical staining of Human Co...