

Product Datasheet

PIK3AP1 antibody (orb30892)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PIK3AP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	90398
Uniprot ID	Q6ZUJ8
NCBI	NP_689522.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

