

Product Datasheet

eEF2K (Phospho-Ser366) Polyclonal Antibody (orb308986)

Description	eEF2K (Phospho-Ser366) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O00418
Expiration Date	12 months from date of receipt.