

Product Datasheet

DARPP-32 (Phospho-Thr34) Polyclonal Antibody (orb309070)

Description	DARPP-32 (Phospho-Thr34) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9UD71
Expiration Date	12 months from date of receipt.