

Product Datasheet

C6orf136 antibody (orb30923)

Description	C6orf136 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Uncharacterized protein C6orf136 protein (1-172AA)
Target	C6orf136
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5SQH8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.