


Product Datasheet

Anti-ITIH2 Antibody (orb309520)

Description	Anti-ITIH2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of fetal br...



Western blot analysis of fetal brain usi...