


Product Datasheet

Anti-JHDM1D Antibody (orb309604)

Description	Anti-JHDM1D Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of fetal br...



Western blot analysis of JHDM1D antibody...