

Product Datasheet

GLT8D1 antibody (Biotin) (orb30970)

Description	GLT8D1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Glycosyltransferase 8 domain-containing protein 1 protein (33-332AA)
Target	GLT8D1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q68CQ7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.