



Product Datasheet

Anti-PLEKHG3 Antibody (orb309705)

Description	Anti-PLEKHG3 Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cell lysat...



Immunohistochemical staining of fetal br...