

Product Datasheet

ATAD3C antibody (FITC) (orb30979)

Description	ATAD3C antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human ATPase family AAA domain-containing protein 3C protein (1-122AA)
Target	ATAD3C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5T2N8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.