

Product Datasheet

Cleaved-Caspase-3 p17 (D175) Polyclonal Antibody (orb310041)

Description	Cleaved-Caspase-3 p17 (D175) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P42574
Expiration Date	12 months from date of receipt.