

Product Datasheet

Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody (orb310045)

Description	Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P55211
Expiration Date	12 months from date of receipt.