

Product Datasheet

NUB1 Polyclonal Antibody (orb310798)

Description	NUB1 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y5A7
Expiration Date	12 months from date of receipt.