

Product Datasheet

KREMEN2 antibody (orb31093)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KREMEN2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48849
Uniprot ID	Q8NCW0
NCBI	NP_001240654.1 , NP_078783.1 , NP_757384.1 , NP_001240655.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

