

Product Datasheet

Ribosomal Protein S9 Polyclonal Antibody (orb310953)

Description	Ribosomal Protein S9 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P46781
Expiration Date	12 months from date of receipt.