

Product Datasheet

HIST1H2AB antibody (orb31102)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Immunogen Synthetic Peptide

Target H2AC4

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 14135

Uniprot ID [P04908](#)


NCBI [NP_066390.1](#), [NP_003504.2](#)

Dilution Range WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of K562, CEM, MDA-...