

Product Datasheet

Tyk 2 Polyclonal Antibody (orb311060)

Description	Tyk 2 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P29597
Expiration Date	12 months from date of receipt.