

Product Datasheet

LCE2A antibody (Biotin) (orb31142)

Description	LCE2A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Late cornified envelope protein 2A protein (1-106AA)
Target	LCE2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5TA79
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.