

Product Datasheet

MEPCE antibody (orb31149)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MEPCE
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74355
Uniprot ID	Q7L2J0
NCBI	NP_001181919.1 , NP_001181920.1 , NP_001181921.1 , NP_062552.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

