



Product Datasheet

Pan Methylated Lysine Monoclonal Antibody (orb311557)

Description	Pan Methylated Lysine Monoclonal Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of HeLa cells usin...