

Product Datasheet

Human Bifunctional 3'-Phosphoadenosine 5'-Phosphosulfate Synthase 1 protein (C-His) [Out of stock] (orb311662)

Description	Recombinant human Bifunctional 3'-Phosphoadenosine 5'-Phosphosulfate Synthase 1 protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl buffer (pH 8.0), 100mM NaCl and 20% glycerol.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Ala624
Uniprot ID	O43252
Expiration Date	6 months from date of receipt.