

Product Datasheet

Human CASB protein (N-His) (orb311678)

Description	Recombinant human CASB protein (N-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl,500mM NaCl,10% Glycerol,1mM DTT,0.1mM PMSF,pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Glu26-Val226
Source	E.coli
Uniprot ID	P05814
Expiration Date	6 months from date of receipt.