

Product Datasheet

Human Complement Factor H protein (C-His) (orb311704)

Description	Recombinant human Complement Factor H protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS,20%Glycerol,5% Trehalose,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 20% Glycerol, 5% Trehalose, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Glu19-Leu449
Source	Human Cells
Uniprot ID	P08603
Expiration Date	6 months from date of receipt.