

## Product Datasheet

### Human FGF-8B protein (Mouse-Human) (orb311760)

|                         |                                                                                                                                                                                                                                                          |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Recombinant human FGF-8B protein<br>(mouse-human)                                                                                                                                                                                                        |
| <b>Endotoxins</b>       | Less than 0.1 ng/ug (1 EU/ug)                                                                                                                                                                                                                            |
| <b>Preservatives</b>    | Lyophilized from a 0.2 µm filtered<br>solution of PBS, pH 7.4.                                                                                                                                                                                           |
| <b>Form/Appearance</b>  | Lyophilized from a 0.2 um filtered<br>solution of PBS, pH 7.4.                                                                                                                                                                                           |
| <b>Storage</b>          | Lyophilized protein should be stored at<br>-20°C, though stable at room<br>temperature for 3 weeks.<br>Reconstituted protein solution can be<br>stored at 4-7°C for 2-7 days. Aliquots<br>of reconstituted samples are stable at -<br>20°C for 3 months. |
| <b>Tag</b>              | tag free                                                                                                                                                                                                                                                 |
| <b>Note</b>             | For research use only                                                                                                                                                                                                                                    |
| <b>Protein Sequence</b> | Gln23-Arg215                                                                                                                                                                                                                                             |
| <b>Source</b>           | E.coli                                                                                                                                                                                                                                                   |
| <b>Uniprot ID</b>       | <b>P55075, P37237</b>                                                                                                                                                                                                                                    |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                                                                                                                                                                           |