

Product Datasheet

Human LIPF protein (C-His) (orb311771)

Description	Recombinant human LIPF protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 25mM TrisHCl, 100mM glycine, 10%glycerol,pH 7.3.
Form/Appearance	Supplied as a 0.2 um filtered solution of 25mM Tris-HCl, 100mM glycine, 10% Glycerol, pH 7.3.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Leu20-Lys398
Source	Human Cells
Uniprot ID	AAI12273.1
Expiration Date	6 months from date of receipt.