

Product Datasheet

Human GMPR protein (C-His) (orb311782)

Description	Recombinant human GMPR protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 40% glycerol, 0.15M NaCl and 1mM DTT, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 40% Glycerol, 0.15M NaCl and 1mM DTT, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Ser345
Source	Human Cells
Uniprot ID	AAH08281.1
Expiration Date	6 months from date of receipt.