

## Product Datasheet

### Human IDE protein (C-His) (orb311815)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Recombinant human IDE protein (C-His)                                                                   |
| <b>Endotoxins</b>        | Less than 0.1 ng/ug (1 EU/ug)                                                                           |
| <b>Preservatives</b>     | Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,0.05%Brij35,10%Glycerol,pH7.5.        |
| <b>Form/Appearance</b>   | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 0.05%Brij35, 10% Glycerol, pH 7.5. |
| <b>Storage</b>           | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.                  |
| <b>Tag</b>               | C-His                                                                                                   |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | 6His tag at the C-terminus.                                                                             |
| <b>Protein Sequence</b>  | Met42-Leu1019                                                                                           |
| <b>Source</b>            | Human Cells                                                                                             |
| <b>Uniprot ID</b>        | <b>P14735</b>                                                                                           |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                          |