

Product Datasheet

Human LCP2 protein (C-His) (orb311882)

Description	Recombinant human LCP2 protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,20% glycerol pH 8.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 200mM NaCl, 20% Glycerol, pH 8.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-T7&C-His
Note	For research use only
Application notes	T7 tag at the N-terminus, 6His tag at the C-terminus.
Protein Sequence	Met1-Pro533
Source	E.coli
Uniprot ID	Q13094
Expiration Date	6 months from date of receipt.