

Product Datasheet

Human Lysozyme g-like protein 2 protein (C-His) (orb311886)

Description	Recombinant human Lysozyme g-like protein 2 protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,10%Glycerol,pH7.2.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.2.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ser20-Phe212
Source	Human Cells
Uniprot ID	Q86SG7
Expiration Date	6 months from date of receipt.