

Product Datasheet

Human MMP-9 protein (C-His) (orb311893)

Description	Recombinant human MMP-9 protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 2mM CaCl ₂ , 150mM NaCl, 0.05% Brij35(w/v), pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 2mM CaCl ₂ , 150mM NaCl, 0.05% Brij35(w/v), pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ala19-Asp707
Source	Human Cells
Uniprot ID	AAH06093.1
Expiration Date	6 months from date of receipt.