

Product Datasheet

Human ACY3 protein (N-His) (orb311903)

Description	Recombinant human ACY3 protein (N-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,100mM NaCl,1mM DTT,10%Glycerol,pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Ser319
Source	E.coli
Uniprot ID	Q96HD9
Expiration Date	6 months from date of receipt.