

Product Datasheet

Human Phosphatidylinositol 5-Phosphate 4-Kinase Type-2 Alpha protein (C-His) (orb311926)

Description	Recombinant human Phosphatidylinositol 5-Phosphate 4-Kinase Type-2 Alpha protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Thr406
Source	Human Cells
Uniprot ID	P48426
Expiration Date	6 months from date of receipt.