

Product Datasheet

Human Esophagus Cancer-Related Gene 2 Protein protein (C-His) (orb312023)

Description	Recombinant human Esophagus Cancer-Related Gene 2 Protein protein (C-His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ser20-Cys85
Source	Human Cells
Uniprot ID	P58062
Expiration Date	6 months from date of receipt.