

Product Datasheet

GABA A Receptor beta 3/GABRB3 Antibody (orb312121)

Description

GABA A Receptor beta 3/GABRB3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human GABRB3 (344-375aa EKTAKAKNDRSKSESNRVDAHGNILLTSLEVH), different from the related mouse and rat sequences by five amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

54116 MW

Uniprot ID


[P28472](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Mouse, Rat, Human

Expiration Date

12 months from date of receipt.



WB analysis of GABRB3 using anti-GABRB3 ...