



Product Datasheet

CD8 antibody (orb312181)


Description	CD8 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse CD8 (51-150/247 aa)
Target	CD8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27 kDa
Uniprot ID	P07725
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:100-500, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lymphoid tissue.



Immunohistochemical staining of mouse lymphoid tissue.



Immunohistochemical staining of mouse spleen tissue.