

Product Datasheet

NAT1 antibody (orb31220)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target NAT1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 33899

Uniprot ID [P18440](#)

NCBI [NP_001153646.1](#), [NP_001153644.1](#),
[NP_001153642.1](#), [NP_000653.3](#), [NP_001153651.1](#),
[NP_001153643.1](#), [NP_001153645.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

