

Product Datasheet

DAP Kinase 1 (phospho-Ser308) antibody (orb312214)

Description	DAP Kinase 1 (phospho-Ser308) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human DAP Kinase 1 around the phosphorylation site of Ser308 (KQ(p-S)VR)
Target	phospho-DAPK1 (Ser308)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	157 kDa
Uniprot ID	P53355
Expiration Date	12 months from date of receipt.