



Product Datasheet

Histone H4 (Acetyl K5) antibody (orb312267)

Description	Histone H4 (Acetyl K5) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthesised acetylpeptide derived from human Histone H4 around the acetylation site of Acetyl Lys5 (RG(Ac-K)GG)
Target	Histone H4 (Acetyl K5)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P62805
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of HeLa cell using Histone ...



Western blot: Lane 1: Human HeLa cell; Lane...