


Product Datasheet

IGJ antibody (orb312289)

Description	IGJ antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human IGJ (31-130/159 aa)
Target	JCHAIN
Preservatives	15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	18kDa
Uniprot ID	P01591
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat sple...