

Product Datasheet

Insulin Receptor Beta (phospho-Tyr1361) antibody (orb312290)

Description

Insulin Receptor Beta (phospho-

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor Beta around the phosphorylation site of Tyr1361 (IP(p-Y)TH)
Target	Phospho-Insulin Receptor Beta (Tyr1361)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68 kDa
CAS Number	3
Uniprot ID	P06213
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.