



Product Datasheet

MAP2 (phospho-Thr1613/1619) antibody (orb312314)

Description	MAP2 (phospho-Thr1613/1619) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MAP2 around the phosphorylation site of Thr1613 and Thr1619 (SR(p-T)PGTPG(p-T)PS)
Target	Phospho-MAP2 (Thr1613/1619)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70/201 kDa
Uniprot ID	P11137
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=2ug/Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of SH-SY5Y. Primary Antibody...



Western blot: Lane 1: SH-SY5Y (Human) Cell...