

Product Datasheet

NFAT4 (phospho-Ser165) antibody (orb312387)

Description	NFAT4 (phospho-Ser165) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NFAT4 around the phosphorylation site of Ser165 (SP(p-S)PA)
Target	phospho-NFATC3 (Ser165)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	116 kDa
Uniprot ID	Q12968
Expiration Date	12 months from date of receipt.