



Product Datasheet

PSMG1 antibody (orb312514)


Description	PSMG1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PSMG1 (21-100/288 aa)
Target	PSMG1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	O95456
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ce...



Western blot analysis of mouse Brain lys...



Immunohistochemical staining of human br...