



Product Datasheet

PPP2R5A antibody (orb312600)

Description	PPP2R5A antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PPP2R5A (301-400/486 aa)
Target	PPP2R5A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	Q15172
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Thymus Lysate.



Western blot analysis of Mouse Ovary Lysate.