

Product Datasheet

PRAGMIN antibody (orb312612)

Description	PRAGMIN antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PRAGMIN/Sgk223 (641-740/1406 aa)
Target	PRAGMIN/Sgk223
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	149 kDa
Uniprot ID	Q86YV5
Expiration Date	12 months from date of receipt.