



Product Datasheet

PADI6 antibody (orb312631)

Description	PADI6 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PADI6 (101-200/694 aa)
Target	PADI6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	77kDa
Uniprot ID	Q6TGC4
Expiration Date	12 months from date of receipt.



Paraformaldehyde-fixed, paraffin embedde...



Paraformaldehyde-fixed, paraffin embedde...