


## Product Datasheet

### PEX6 antibody (orb312746)

|                            |                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------|
| <b>Description</b>         | PEX6 antibody                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                |
| <b>Reactivity</b>          | Mouse                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                          |
| <b>Tested Applications</b> | IHC-P                                                                                 |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human PEX6 (751-850/980 aa)             |
| <b>Target</b>              | PEX6                                                                                  |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                 |
| <b>Form/Appearance</b>     | Liquid: Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| <b>Concentration</b>       | 1mg/ml                                                                                |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.      |
| <b>Note</b>                | For research use only                                                                 |
| <b>Isotype</b>             | IgG                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                            |
| <b>Purity</b>              | affinity purified by Protein A                                                        |
| <b>MW</b>                  | 104kDa                                                                                |
| <b>Uniprot ID</b>          | <b>Q13608</b>                                                                         |
| <b>Dilution Range</b>      | IHC-P=1:100-500                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                       |



IHC-P image of rat brain tissue using PE...